

TECHNICAL DATA SHEET

PRIMUS 14" / 16" Natural PF

STERILE LATEX SURGICAL GLOVES BEADED 14" / 16" (350 – 400 mm)

POWDER FREE CHLORINATED / POLYMER COATED

Physical Dimensions

Size	Length (mm) Min	Palm Width	Thickness (mm) Minimum		
			Cuff	Palm	Finger
6.0		76 ± 6			
6.5		83 ± 6			
7.0		89 ± 6			
7.5	350	95 ± 6	0.13	0.19	0.21
8.0		102 ± 6			
8.5		108 ± 6			

Physical Properties

Characteristics	Before Ageing	After Ageing
Tensile Strength (Mpa) Min	24	18
Ultimate Elongation (%) Min	750	560
Stress at 500% elongation (Mpa) Max	5.5	NA
Force at Break (N) Min	12	9

Protein Content <50 µg/dm²

Formulation

Curing Agent	Sulphur
Accelerators	Zinc diethyl dithio carbamate, Zinc mercapto benzthiazole
Activators	Zinc Oxide
Antioxidants	Wingstay L (Condensation product of p-cresol & dicyclopentadiene
Stabiliser	Potassium Hydroxide
Mould Cleaner	Nitric Acid
Coagulant	Calcium Nitrate

Performance Requirements

Characteristics	Related Defects	Inspection Level	AQL
Water tightness (Final Release)	Holes	G-1 Single Sampling Plan	1.5

Primary Material	Latex
Colour	Pale Yellow
Shape	Anatomic
Cuff	Beaded
If Sterilised	Ethylene Oxide / GAMMA
Certificates	ISO 13485 : 2003 , CE Mark
Complies to	EN 455 , ASTM 3577



Packing	Size	Product Code
1 Pair in a Paper Wallet, then to a paper pouch	6.0	1746226 / T
50 Pouches (Pairs) in 1 Inner Box	6.5	1746226.5 / T
8 Inner Boxes (400 Pairs in 1 Master Shipper Carton)	7.0	1746227 / T
	7.5	1746227.5 / T
Pouch Size : 15 x 30 Cm Peel Open Type	8.0	1746228 / T
	8.5	1746228.5 / T
	9.0	1746229 / T

Manufactured By : PRIMUS Gloves Pvt Limited, Plot No 14-A, Cochin Special Economic Zone, Kakkannad, Cochin 682037, Kerala, INDIA
 Tel : + 91 484 2413063, 2413076 ; Fax : + 91 484 2413089
 Email: .marketing@primusgloves.com; primusgloves@gmail.com
 Web : www.primusgloves.com



Issue date : 01.09.2010