



PEEL DOWN TO OPEN

PROTAC PF

Powder Free

Latex

Powder Free Sterile Surgical Gloves

CAUTION: This product contains natural rubber latex which may cause allergic reactions.
Sterile if pack undamaged or unopened.

- The latest generation in powder free surgical gloves
- New synthetic inner coating for :-
 - a. Easy donning even with damp hands
 - b. Easy intra-operative glove changes.
 - c. Reduced contact with latex allergens & Residual process chemicals.
 - d. New anatomical shape & excellent grip
 - e. Comfort level equivalent to powdered gloves (no hand fatigue during long operation)
 - f. Reduced latex protein content & residual process chemicals
 - g. Beaded cuff ensure a secure hold on the surgical gown



Single use only



1 Pair

STERILE EO

Sterilized by Ethylene Oxide



US FDA
510 K




Primus Gloves Pvt Limited

Plot No. 14-A, Cochin Special Economic Zone, Kakkanad, 682037, Kochi, Kerala.


www.primusgloves.com

Customer Care No.: 0484-2413076, sterimedmktg@primusgloves.com, sterimedcochin@gmail.com

Made in India

Mfg.Date  :

UseBefore  :

BatchNo.  :

Size :

EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

PRIMUS
Your Partner In Protection