



Cut to open

Protac

7

Powdered

Latex

STERILE

Disposable Surgical Rubber Gloves

CAUTION: This product contains natural rubber latex which may cause allergic reactions. Sterile if pack undamaged or unopened.

Warning: Surface powder shall be removed aseptically prior to undertaking operative procedures in order to minimize the risk of adverse tissue reactions.

- Micro rough surface for secure, positive grip.
- Natural rubber latex for excellent fit, feel & comfort
- Lightly powdered for ease of donning.
- AQL 1.5 - global standard for low pinhole incidence.
- Beaded Cuff for Excellent gown grip



Single use only



1 Pair

STERILE EO

Sterilized by Ethylene Oxide



US FDA
510 K



PRIMUS GLOVES PVT LIMITED.
14-A, CSEZ, Cochin - 682037

Customer Care No : 0484-2413076
sterimedmktg@primusgloves.com, sterimedcochin@gmail.com

Mfg. Date  :

Use Before  :

Batch No.  :

Size : 7

Made in India

EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

PRIMUS

Your Partner In Protection

