



Surgical Gloves Powder Free

PRIMUS Gloves Pvt Limited

Plot No. 14A, CSEZ, Kakkanad, Cochin 682 037, Kerala, India

Tel: +91 484 2413063, 2413076 Fax: +91 484 2413089

Email: marketing@primusgloves.com Website: www.primusgloves.com

Surgical Gloves Powder Free

Natural Rubber Latex, Micro Textured, Anatomical in Shape



PRIMUS Surgical Gloves Powder FREE is the result of vast experience and knowledge of the materials. Our quality and production team has analyzed the surgical gloves available in the market and collected feedback from diverse customers and designed a curved former adapting to the hand anatomy ensuring perfect fit just like a second skin, preventing fatigue during long operating procedures. Our formers have micro rough surface to ensure excellent grip. The whole manufacturing process is carried out under stringent process using state of the art automatic glove dipping machine and quality control to ensure a good glove with excellent skin tolerance. We have Chlorinated Powder FREE Surgical Glove and Polymer Coated Powder FREE Surgical Gloves in our range. The gloves are sterilized either by ETO (Ethylene Oxide Gas) or by Gamma Radiation as per customer requirement. Packaging is based on specifications requested by the customer.

- Manufactured from High quality Natural Rubber Latex
- In-process control covering all aspects of manufacturing ensure a consistent high quality product.
- Acceptable Quality Level (AQL) of 1.5. Single Sampling Plan; level G1 as per ISO 2859-1
- In-line leaching process reduce content of soluble proteins
- Avoidance of thiurams for lower allergenic potential
- Manufactured and packed under clean and hygienic environment ensuring Low Bio-Burden on the gloves
- Skin-friendly presentation
- Powder FREE Gloves are manufactured in sophisticated automatic Chlorinator machines for consistent quality
- Consistently low protein levels of 50 µg / dm² or less
- Good elongation and full anatomic shape provide little hand fatigue and lasting comfort in use even during long procedures
- Micro rough surface for that extra grip
- Beaded cuff for secure fit
- Length 280-285 mm
- Glove Size - 6, 6½, 7, 7½, 8, 8½ & 9

Usage:

- Medical / Hospital Use
- Laboratories/Pharmaceutical Industries.

Packing:

Pouch Size	Pouch Type	Packing	Container Quantities	
15 x 30 cm	Peel Down to open	1 pair gloves in a paper wallet, then to a paper pouch; 50 Pairs in 1 Inner Box, 8 Inner Boxes (200 Pairs) in 1 Master Shipper Carton	20 FCL	160,000 Pairs
			40 FCL	336,000 Pairs

Size	Product Code
6.0	1244326/T
6.5	1244326.5/T
7.0	1244327/T
7.5	1244327.5/T
8.0	1244328/T
8.5	1244328.5/T
9.0	1244329/T

PRIMUS Gloves Pvt Limited

Plot No. 14A, CSEZ, Kakkanad, Cochin 682 037, Kerala, India

Tel: +91 484 2413063, 2413076 Fax: +91 484 2413089

Email: marketing@primusgloves.com Website: www.primusgloves.com